


# Medical Alarm Signal Analyzer

## IEC 60601-1-8 Auditory and Visual Alarm Testing

Model: KP-108

### 1 Product Overview

The KP-108 Medical Alarm Signal Analyzer is designed for evaluating auditory alarm signals and visual alarm indicators of medical electrical equipment. It supports IEC 60601-1-8 related alarm testing with waveform acquisition, parameter analysis, and automatic report generation.



### Key Features

- ✓ Auditory alarm signal analysis
- ✓ Visual alarm indicator analysis
- ✓ Time-domain waveform and FFT spectrum
- ✓ Automatic PASS/FAIL report output
- ✓ USB and LAN data export

### Typical Applications

- Patient monitors
- Infusion pumps
- Ventilators
- Anesthesia systems
- Dialysis systems

### Standards and Reports

#### Standard Reference

IEC 60601-1-8 related testing

#### Test Types

Auditory and visual alarm signals

#### Report Output

Waveform, parameters, summary

# Testing Capabilities and Report Output

KP-108 Medical Alarm Signal Analyzer

## 1 What the Analyzer Measures

### Auditory Alarm Signals

- Time-domain waveform
- FFT spectrum
- Pulse duration and interval
- Dominant frequency and harmonics
- Priority recognition support

### Visual Alarm Indicators

- Flashing frequency
- Duty cycle
- Rise and fall time
- Optical response level
- Wavelength range

## 2 Typical Test Workflow



### 3. Key Benefits

- ✔ Supports IEC 60601-1-8 related evaluation
- ✔ Improves test repeatability and efficiency
- ✔ Clear report output for laboratory documentation
- ✔ Suitable for R&D, quality control, and compliance testing


### 4. Typical Deliverables

- Test waveform screenshots
- Parameter tables
- PASS/FAIL summary
- Exportable report records

# Auditory Alarm Signal Analysis Report


Time-domain waveform, FFT spectrum, alarm parameters and PASS/FAIL summary

## 1. Time Domain Waveform



Input Range:  $\pm 2.0$  Pa  
 Sampling Rate: 48.0 kHz  
 Record Length: 5.000 s  
 RMS Level: 0.281 Pa  
 Peak Level: 0.932 Pa  
 Crest Factor: 3.31  
 DC Offset: -0.001 Pa

## 2. Frequency Domain (FFT Spectrum)



FFT Size: 16384  
 Window: Hann  
 Frequency Res.: 2.93 Hz  
 A-Weighting: On  
 Ref. Level: 1.000 Pa

## 3. Alarm Signal Parameters

Parameter	Value	Unit
Pulse Duration ( $T_{on}$ )	220	ms
Pulse Interval ( $T_{off}$ )	780	ms
Pulse Period ( $T_{total}$ )	1000	ms
Rise Time (10-90%)	12.4	ms
Fall Time (90-10%)	18.7	ms
Burst Count (in record)	5	bursts
Dominant Frequency ( $F_{1,1}$ )	1520	Hz
Harmonic Content (2nd-5th)	-13.2	dB (avg.)
Overall Level (A-weighted RMS)	0.281	Pa
Priority Recognition Result	<b>High Priority</b>	(recognized)

## 4. Analysis Summary

**Test Result:** **PASS**  
**Standard/Template:** IEC 60601-1-8:2006 (High Priority)  
**Compliance:** Within Limits  
**Analysis Time:** 2025-05-20 14:32:18


### Analysis Notes:

- Pulse pattern: 5 bursts, regular.
- Measured pulse duration and interval are within the typical limits for High Priority alarms.
- Dominant frequency at 1520 Hz with harmonics up to 5th order.
- Signal level stable; no clipping or distortion detected.
- Priority recognition algorithm result: High Priority (confidence 96%).


# Visual Alarm Indicator Analysis Report

Flashing waveform, pulse train, visual alarm parameters and status summary

## 1. Flashing Signal Waveform (Optical Response vs. Time)



## 2. Pulse Train View (Visual Indicator)



## 3. Visual Alarm Indicator Parameters

Parameter	Measured Value	Specification	Result
Flashing Frequency	1.00 Hz (60.0 fpm)	0.5 – 2.0 Hz (30 – 120 fpm)	PASS
Duty Cycle	49.8 %	30 – 70 %	PASS
Rise Time (10% – 90%)	12.3 ms	≤ 100 ms	PASS
Fall Time (90% – 10%)	11.7 ms	≤ 100 ms	PASS
Wavelength Range (Peak)	629 nm	620 – 635 nm	PASS
Optical Response (Peak)	101.2 mW/m <sup>2</sup>	≥ 20 mW/m <sup>2</sup>	PASS
Optical Response (Min)	0.18 mW/m <sup>2</sup>	≤ 5 mW/m <sup>2</sup>	PASS

All measurements performed in accordance with IEC 60601-1-8:2020, Clause 8.6 (Visual Alarms).

## 4. Status Summary

Overall Test Result: **PASS**

### Analysis Notes:

- Visual alarm indicator flashes within specified frequency range.
- Duty cycle is within allowable limits.
- Rise and fall times meet response requirements.
- Peak wavelength and optical response are within specification.
- No irregularities detected in flashing waveform.

Test Mode:	Continuous Capture	Detector:	Si Photodiode
Measurement Range:	0 – 200 mW/m <sup>2</sup>	Optical Probe S/N:	PD-1000-2A-0187
Ambient Light:	12.5 lux	Calibration Date:	Mar 12, 2025
Test Distance:	1.0 m	Next Calibration Due:	Mar 12, 2026